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

### Delta Scaling:

Scale setpoint based on which joint is going to take the longest to achieve. Allows all joints to arrive simultaneously.

### Velocity Scaling:

'Speed Ratio' describes the overall velocity scaling.

### High Speed Scaling:

High speed scaling reduces execution speed when commanded speed exceeds a high speed velocity threshold **and** the arm's high-speed collision links are in collision.

### Collision Avoidance:

Applies offsets to joint commands based on depth of intersection between arm collision geometries and the opposing arm or torso.


### Joint Limit Clipping:

If the joint command is beyond limits, clip the command to respect joint limits.

### Collision Detection:

If collision (impact) is detected, set position command to hold current compensating for the impact.

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**High Speed Scaling:**

High speed scaling reduces execution speed when commanded speed exceeds a high speed velocity threshold **and** the arm's high-speed collision links are in collision.

**Collision Avoidance:**

Applies offsets to joint commands based on depth of intersection between arm collision geometries and the opposing arm or torso.

**Joint Velocity Clipping:**

Limits joint velocity command to not exceed maximum joint velocities.


**Joint Limit Check:**

Validates that resulting joint position will be within joint limits. If not, no velocity will be commanded to any joint.

**Collision Detection:**

If collision (impact) is detected, set position command to hold current compensating for the impact.

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**Gravity/Spring Compensation:**

The joint torque command is applied in addition to the gravity and S1 spring compensation torques.

**Joint Torque Scaling:**

Scales all joint torques if a torque command exceeds the maximum allowable torque for that joint. This scaling ratio is defined as  $\text{torque\_max} / \text{torque\_command}$ .